Rollator

Your Home Rehabilitation Partner

Z21

V4.1

@VOCiCMedical

@VOCiCMedical

@VOCiCMedical

@VOCiCMedicalOfficial

EN CAUTION

- This product is designed for its intended purpose only. Please read this manual thoroughly before first use. It will guide you on how to assemble and operate the Z21 Rollator.
- Before assembly and use, ensure that all components and accessories are complete and undamaged.
- If you find any issues, do not use the product and contact us immediately at support@vocic.com.
- Please follow the instructions carefully during assembly. Do not attempt installation without consulting the instructions, as improper assembly may damage components and create safety hazards.
- Before use, assess whether the user is suitable for this product based on factors such as weight, height, and mobility.
- The product should not be used by individuals who are not suitable, as this may lead to injury. Disabled users should use the product only when accompanied by a companion or with an emergency alert device.
- Do not disassemble or modify any parts of the product. Unauthorized modifications may damage the product and present safety risks.
- This package contains a plastic bag. To prevent suffocation, keep the plastic bag out of reach of children. The plastic bag is not a toy, and should never be given to infants or toddlers.
- This product is not a toy. Do not allow children to operate or use it.

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Manufacturer: Zhejiang Mateside Medical Devices Technology Co., Ltd.
Address: Room 501, Unit C, Building 1, 98 Luze Jie, Cangqian Street, Yuhang District, Hangzhou, Zhejiang, CHINA
support@vocic.com

EC REP: Name: **SUNGO Europe B.V.** +31 (0) 10 3034500 ec.rep@sungoglobal.com
Address: Fascinatio Boulevard 522, Unit 17, 2909VA Capelle aan den IJssel, The Netherlands

UKRP: Name: **SUNGO Certification Company Limited** +44- (0) 7868 234015 ukrp@sungoglobal.com
Address: 3rd floor, 70 Gracechurch Street, London, EC3V 0HR

**MADE IN CHINA**
Corresponding to the UDI label

Expanded Width: 59cm / 23.2in

Do not move rollator while seated.
For Indoor & Outdoor Use

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











Please lock the rollator before sitting!

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
IV. Customer Support & After-Sales Service

Safety Instructions

⚠ Important Notice

- Thank you for choosing our product. Throughout this guide, we will use the model name Z21 as a shortened reference for the full product name, Z21 Rollator. Please keep this in mind as you review the information provided.
- For an optimal experience, please read this manual carefully before using the Z21 for the first time. It will guide you through the assembly and proper use of the Z21.
- Please retain this manual and make it available to anyone who may use or come into contact with the product for reference.
- This product is intended solely for its designed purpose. Do not use it as a wheelchair or for any other similar purposes.

⚠ Operation Warnings

- Before use, confirm the maximum weight capacity of this product is **158 kg (350 lb)**. Do not exceed this limit to prevent damage to the product or potential safety hazards.
- Prior to use, ensure all connectors and fasteners are securely installed to avoid product damage or safety risks caused by improper installation.
- Test the braking system before use: When the handbrake is engaged, the rear wheels should be locked and unable to move. When the handbrake is released, the wheels should rotate freely.
- Before each use, check the tightness of screws and knobs, and retighten them if necessary.
- During use, folding, or unfolding, keep hands and objects away from movable parts such as the area beneath the seat, the gaps between, or leg tubes to prevent injury or other safety risks.
- After folding and storing the product, ensure it is fully and correctly unfolded before reuse to prevent potential safety risks.
- The seat provided with this product is suitable for sitting and resting. However, always engage the handbrake to ensure the rear wheels are locked before using the seat.
- Do not move the product while someone is seated. This product is not designed to function as a transportation device; improper use may result in injury or accidents.
- Avoid using this product on slippery surfaces and do not use the seat function on steep slopes or uneven terrain.

⚠ General Safety Warnings

- Before installation and use, please ensure that all components and accessories are complete, intact, and free from damage.
- If there are any issues, do not use the product and contact us promptly at support@vocic.com.
- When assembling the product, strictly follow the installation steps provided in the manual. Do not attempt self-assembly to avoid potential damage to the components or safety risks caused by improper installation.
- Before use, evaluate whether the user is suitable for this product, considering factors such as weight, height, and mobility. Individuals who are not suitable should refrain from using the product to avoid injury. Persons with disabilities should only use the product with a companion present or in environments equipped with alarm devices. Do not disassemble or modify any parts of this product without authorization, as this may cause damage and pose safety risks.
- This package contains plastic bags. To avoid the risk of suffocation, keep plastic bags out of reach of children. Plastic bags are not toys, do not allow infants or young children to handle them.
- This product is not a toy and should not be operated or used by children.

⚠ Maintenance and Care

- Users are advised to regularly inspect the device to ensure it is functioning properly.
- During use or maintenance, do not clean this product with disinfectants or other cleaning agents. If necessary, use a soft cloth or brush with a small amount of water for cleaning. After cleaning, wipe the product dry to ensure it remains in a dry state.
- When cleaning the product, do not spray it directly with a showerhead or similar device. This may wash away the lubricant from moving parts, reducing the product's lifespan.
- For long-term storage, ensure the product is completely dry before packaging to prevent rust or other issues caused by prolonged storage.

🔔 Notice

If you notice any visible damage, unusual sounds, or malfunctions, please stop using the product immediately. Refer to the **Troubleshooting** section for assistance. If the issue persists, please contact us immediately at support@vocic.com. Thank you.

Product Information

Product Name: Z21 Rollator

Purpose

Specifically designed to provide stable support for individuals with disabilities, the elderly, or those with weakened lower limbs. Assists with walking and standing, helping users perform daily activities more safely and independently. Suitable for various settings including home use and outdoor environments.

Technical Specifications

Expanded Width:	59 cm / 23.2 in
Expanded Length:	65-70 cm / 25.6-27.5 in
Folded Dimensions:	35×59×80-98cm / 13.8×23.2×31.5-38.6in
Handle Height Range:	80-93 cm / 31.4-36.6 in
Seat Dimensions:	36×31 cm / 14.2 × 12.2 in
Seat Height:	52-56 cm / 20.4 × 22 in
Front Wheel Size:	20 cm / 8 in
Rear Wheel Size:	20 cm / 8 in
Net Weight:	9.8 kg / 21.6 lb
Weight Capacity:	158 kg / 350 lb
Frame Material:	Iron
Wheel/Handle Material:	PVC
Seat Material:	PU
Backrest Material:	Oxford Fabric+PE
Storage Bag Material:	Oxford Fabric

Troubleshooting

Fault Description	Fault Check	Solution
The rear wheels continue to rotate even when the handbrake is pulled or pressed.	The brake may be too loose.	Please refer to Part IX of the manual for brake tension adjustment instructions. If this solution does not resolve the issue, feel free to contact us at support@vocic.com .
Difficulty pushing the device or unusual noises during operation.	The wheels may be installed incorrectly, left and right reversed.	Check if the wheels are installed correctly, following the instructions in the manual. If the installation seems correct, please contact us at support@vocic.com for assistance.
	The knobs or screws may be installed incorrectly.	
	The parts may not be properly tightened.	Check if all parts are properly tightened. If the installation appears correct, please contact us at support@vocic.com for assistance.

After-Sales Service

After-Sales Contact

• Standard After-Sales

If you encounter any of the following issues, please feel free to contact us at support@vocic.com. To ensure optimal service, we will provide an effective solution within 24 hours. Thank you for your support and understanding!

• Scope of After-Sales Service

			
Missing Parts	Damaged Components	Product Malfunction	Usage Guidance Service

Return and Exchange Policy

If you need to return or exchange your product, please feel free to contact us at support@vocic.com. In support of environmental protection and to reduce unnecessary carbon emissions and resource waste, we kindly ask that you retain and use the original packaging for returns and exchanges. We will provide an effective solution within 24 hours. Thank you for your support and understanding!

To qualify for a return or exchange, please ensure the following conditions are met:

- The product remains in its original packaging and in good condition, without any damage or contamination.
- All accessories, manuals, and other included items are complete and intact.
- A valid proof of purchase is provided.

Note: If the product does not meet these return and exchange conditions, the request will be declined, and any related costs will be the responsibility of the buyer.

Warranty Information

• Warranty Period

This product is covered by a 3-year warranty starting from the date of purchase. If any faults arise during the warranty period due to material or craftsmanship issues, we will provide free replacement services.

• Warranty Coverage

- The warranty covers the main functional components of the product.
- The warranty is valid only for product quality issues arising from normal use, and does not cover damage caused by human error.

How to Contact Us

Reason for Contact:

☐ After-Sales Support

☐ Return and Exchange

☐ Warranty

Required Information:

Product Name / Model number /

Purchase Platform:

Order Number:

Email:

Issue Description:

Email: support@vocic.com

Disclaimer

Thank you for choosing our product! To ensure safe and proper use, please read the following disclaimer carefully: Before using this product, you should thoroughly understand the user instructions and operation methods provided. Any damage or injury caused by improper use or failure to follow the instructions will be your responsibility. By using this product, you acknowledge and agree to assume all risks associated with its operation.

Under normal use and appropriate maintenance, the recommended service life of this device is three years (note that the recommended service life does not constitute a warranty period. The warranty covers defects under normal usage conditions, while the recommended service life represents the expected optimal duration of use for the device).

Our company will not assume responsibility for any damages caused by events including but not limited to war, governmental policies, fire, earthquakes, floods, typhoons, tsunamis, power outages, or strikes. We shall also not be liable for any anticipated benefits, reliance benefits, indirect damages, or special or incidental events (including, but not limited to, economic losses, data loss, or personal injury) resulting from the use of this product.

This product is provided as is, without any explicit or implied warranties. Our company does not guarantee that this product will be suitable for all usage scenarios or meet specific user requirements. Please assess this product's suitability based on your own needs before purchase and use. The warranty services we provide are limited to quality issues inherent to the product itself. You assume all risks and responsibilities during use. For detailed terms of the warranty, please refer to the Warranty Guide.

Before using the product, carefully inspect its condition to ensure it is in good working order. If any abnormalities are detected during use, please immediately stop using the product and contact our company or an authorized service provider.

Our company reserves the right to modify, update, and provide final interpretations of this disclaimer. If you have any questions or concerns, please contact customer service for further assistance.

This disclaimer applies to all users who purchase and use our products. By using this product, you acknowledge and accept the terms of the above disclaimer.

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